

CYIENT

Deep Learning-Based Testicular

Ultrasound Quantification

Automated segmentation delivers faster, more reproducible
scrotal ultrasound assessment for infertility evaluation.



CONTENTS

Abstract	03
Introduction	03
Clinical Workflow Challenges	03
Technology Overview	04
Automated 2D Quantification Parameters	04
Workflow Benefits	05
Technical Considerations	05
Intended Use and Limitations	05
The Cyient Thought Board	06
Conclusion	07
About the Authors	07
About Cyient	08

Abstract

Manual testicular measurement using electronic calipers is standard in scrotal ultrasound, but it can be time-consuming and operator-dependent, with variation arising from image plane selection, caliper placement, and image quality. This whitepaper presents a deep learning-based approach to automated 2D testicular quantification that segments the testicular boundary and derives objective morphometric parameters, including area, major axis, minor axis, equivalent diameter, and perimeter. By standardizing measurement across cases, the technology helps reduce inter-operator variability, speed up reporting, and improve reproducibility for clinical follow-up and research – as a decision-support tool that complements, rather than replaces, clinical interpretation.


Role of Ultrasound in Male Infertility Assessment

Scrotal ultrasound is routinely used to assess testicular size, morphology, echotexture, and related clinical findings in male infertility evaluation. In standard practice, testicular measurements are performed manually by placing electronic calipers along the longest visible testicular axis and the orthogonal short axis. While widely accepted, this approach can be time-consuming and operator-dependent, particularly when multiple planes, bilateral testes, serial follow-up scans, or research datasets are involved.

Manual testicular measurement can be time-consuming and operator-dependent.

Limitations of Manual Testicular Measurement

Manual caliper measurement can vary depending on:



Even small differences in caliper placement may produce measurable variation, particularly in borderline cases or longitudinal follow-up, and repeated manual measurement increases reporting time in high-throughput clinics and research studies.

Deep Learning-Based Testicular Segmentation

Deep learning-based automated 2D testicular quantification uses a trained segmentation model to automatically identify the testicular boundary on 2D grayscale ultrasound images. A measurement engine then computes quantitative morphometric parameters from the generated mask. The workflow follows four steps:

Automated segmentation reduces measurement to a single image review step.





Figure 1. Deep learning-based testicular segmentation and measurement interface.

Automated 2D Quantification Parameters

The measurement engine calculates area (mm² or cm²), the longest Feret or fitted-axis diameter (major axis), the orthogonal short-axis measurement (minor axis), boundary length (L), and shape descriptors such as eccentricity/aspect ratio, using image pixel spacing to report values in physical units. These automated measurements provide a consistent quantitative representation of the testis on the selected 2D plane and, when integrated into standard imaging workflow, can reduce the need for repeated manual caliper placement across bilateral, serial, and multi-image studies.

Operational Benefits in Ultrasound Measurement





Figure 2. Comparison of manual and automated testicular measurements.

THE CYIENT THOUGHT BOARD

Q: Why is manual testicular measurement inconsistent?

Depends on caliper placement and operator experience

Q: What does the AI system do?

Automatically segments the testis and measures it

Q: What parameters does it measure?

Area, major axis, minor axis, perimeter, eccentricity

Q: How does it improve workflow?

Cuts measurement time and reporting variability

Q: Can clinicians still review the results?

Yes – operators confirm or correct segmentation

Q: Does it replace clinical judgment?

No – physician interpretation remains essential

Conclusion

Deep learning-based 2D ultrasound testicular quantification enables rapid, standardized, and reproducible measurement of testicular morphology. By automatically segmenting the testis and extracting area, major axis, minor axis, and related parameters, the technology can reduce manual measurement time and minimize operator-dependent variability. It supports efficient scrotal ultrasound workflow, improves consistency in measurement reporting, and provides objective imaging biomarkers for clinical review, research, and longitudinal assessment – always as a decision-support and workflow-assistance tool, with final interpretation remaining under qualified clinical supervision.

Automated
measurement
supports –
but does not
replace –
clinical
interpretation.

About the Authors



Srinivas Rao Kudavelly is Senior Director in the Technology Office at Cyient Limited, Bangalore, Karnataka, where he leads the development of AI and deep learning-based medical imaging solutions.



Dr Rama Raju is Director at Krishna IVF Clinic, Vizag, Andhra Pradesh, bringing clinical expertise in reproductive medicine and male infertility assessment.

About Cyient

Cyient is a global lifecycle engineering company powering mission-critical industries from design to aftermarket across products, plants, and networks. With 300+ customers and 15,000+ associates across 30+ countries, we combine human expertise, deep domain knowledge, and AI to deliver intelligent engineering solutions across Aerospace and Rail, Automotive and Mobility, Connectivity, Energy, Healthcare, Mining, Utilities, and beyond.

To learn more, visit www.cyient.com



Contact Us

North America Headquarters

Cyient, Inc.
99 East River Drive
5th Floor
East Hartford, CT 06108
USA
T: +1 860 528 5430
F: +1 860 528 5873

Europe, Middle East, and Africa Headquarters

Cyient Europe Limited
Apex, Forbury Road,
Reading
RG1 1AX
UK
T: +44 118 3043720

Asia Pacific Headquarters

Cyient Limited
L14 3 Parramatta Square,
153 Macquarie St,
Parramatta NSW 2150,
Australia
T: +61 2 8887 8600

Global Headquarters

Cyient Limited
Plot No. 11
Software Units Layout
Infocity, Madhapur
Hyderabad - 500081
India
T: +91 40 6764 1000
F: +91 40 2311 0352

Follow us on:  